

PMA Hybrid Series



Hybrid Photomultiplier Detector Assembly

- Detection efficiency up to 45 % at 500 nm (cathode dependent)
- Instrument response down to 50 ps (FWHM, cathode dependent)
- Negligible afterpulsing
- Internal HV power supply and pre-amplifier
- Active temperature stabilization
- Shutter and overload protection
- Active sensor area: 3 mm, 5 mm or 6 mm (cathode dependent)
- Analog output provides positive voltage proportional to count rate



Applications

- Fluorescence lifetime measurements
- Fluorescence Lifetime Imaging (FLIM)
- Fluorescence Correlation Spectroscopy (FCS, FCCS, FLCS, FLCCS, 2-focus FCS)
- Ranging (LIDAR, OTDR)
- Laser scanning microscopy (confocal / two-photon)

The PMA Hybrid is a compact, single photon sensitive detector based on a fast hybrid photomultiplier tube (R10467 from Hamamatsu) with Peltier cooler to reduce the dark count rate. It features two signal outputs: a negative voltage pulse for timing and counting applications and an analog positive output voltage proportional to the count rate that can be connected to e.g. A/D converters. The detector includes a high voltage power supply, a pre-amplifier with overload protection, and an emergency shut down procedure if the detector count rate reaches a critical limit. Overload protection, high voltage set-up, and temperature regulation are calibrated at PicoQuant and do not require any user adjustment. The detector's CAN interface is compatible with all PicoQuant systems.


Five different photocathodes can be incorporated into the PMA Hybrid to meet the user's needs. The PMA Hybrid 06 and 07 are sensitive in the UV and in the blue spectral region. High detection efficiencies in the visible spectrum (up to 45 % at 500 nm) are achieved with the PMA Hybrid 40. For applications using light at longer wavelengths of the visible spectrum and in the NIR, the PMA Hybrid 42 and 50 are the best choices.

The PMA Hybrid is encased in a nickel coated aluminum housing to achieve high level of RF shielding and protection against interference from other devices. The built-in pre-amplifier is specially optimized for timing sensitive applications such as Time-Correlated Single Photon Counting (TCSPC).


All PMA Hybrids have very good timing resolution that can even reach values down to 50 ps (FWHM) for the blue sensitive version. In contrast to other detector types, afterpulsing is negligible. This special feature makes the PMA Hybrid especially suited for e.g., Fluorescence Correlation Spectroscopy (FCS), where the afterpulsing peak at early lag times often complicates the analysis of the autocorrelation function.

The PMA Hybrid interfaces directly with all PicoQuant TCSPC products such as the PicoHarp 300 or HydraHarp 400. Due to its large active area, the detector can be connected to spectrometers like the FluoTime 200 or FluoTime 300 from PicoQuant. It can also be attached to Laser Scanning Microscopes in Non-Descanned Detection (NDD) set-ups via the C-mount adapter. Integration in descanned detection mode or other systems, e.g., the confocal time-resolved microscope MicroTime 200 from PicoQuant is also possible. The output signal of the PMA Hybrid is accessible through a standard 50 Ohms SMA connector. The module only needs a 12 V DC supply line, which is included in each PMA Hybrid delivery.


Spectral response



Measurement examples



FCS curve (autocorrelation) of a 1 nm ATTO 488 solution. No afterpulsing peak visible.



Fluorescence decay of a Anthracen solution in EtOH shows clean response of the detector.

Specifications

Cathode type	-06	-07	-40
Electrical parameters			
Wavelength range	220 to 650 nm	220 to 850 nm	300 to 720 nm
Dark counts (cooled, typ. value)	< 100 cps	< 200 cps	< 700 cps
Transit time spread (FWHM, typ. value)	< 50 ps	< 50 ps	< 120 ps
Recommended max. count rate	10 MHz		
Overload shutdown at	80 MHz (with cw excitation, lower values at different conditions)		
Single electron response width (typ. value)	600 ps		
Pulse rise / fall time (typ. value)	400 ps		
Signal output (Timing)			
Connector	SMA female		
Impedance	50 Ohms		
Polarity	negative		
Signal output (Analog)			
Connector	SMA female		
Impedance	>1k Ohms		
Polarity	positive		
Max. Out	+10V (corresponds to 50 Mcps)		
Time constant of the amplifier	20 µs		
Power supply			
Input	12 V DC		
Dimensions			
Detector area diameter	6 mm	6 mm	3 mm
Housing (w × d × h)	80 × 170 × 60 mm		
Optical adapters	C-mount, mounting holes for FluoTime 200 / FluoTime 300 spectrometers		

Cathode type	-40 mod	-42	-50
Electrical parameters			
Wavelength range	300 to 720 nm	300 to 870 nm	380 to 890 nm
Dark counts (cooled, typ. value)	< 4000 cps	< 500 cps	< 1000 cps
Transit time spread (FWHM, typ. value)	< 120 ps	< 130 ps	< 160 ps
Recommended max. count rate	10 MHz		
Overload shutdown at	80 MHz (with cw excitation, lower values at different conditions)		
Single electron response width (typ. value)	600 ps		
Pulse rise / fall time (typ. value)	400 ps		
Signal output (Timing)			
Connector	SMA female		
Impedance	50 Ohms		
Polarity	negative		
Signal output (Analog)			
Connector	SMA female		
Impedance	>1k Ohms		
Polarity	positive		
Max. Out	+10V (corresponds to 50 Mcps)		
Time constant of the amplifier	20 µs		
Power supply			
Input	12 V DC		
Dimensions			
Detector area diameter	5 mm	3 mm	3 mm
Housing (w × d × h)	80 × 170 × 60 mm		
Optical adapters	C-mount, mounting holes for FluoTime 200 / FluoTime 300 spectrometers		

Please check our website for updated information.

All Information given here is reliable to our best knowledge. However, no responsibility is assumed for possible inaccuracies or omissions. Specifications and external appearances are subject to change without notice. Trademarks or corporate names are used for explanation and identification, to the owner's benefit and without intent to infringe.

PicoQuant GmbH, May 2017



PicoQuant GmbH
Rudower Chaussee 29 (IGZ)
12489 Berlin
Germany

Phone +49-(0)30-1208820-0
Telefax +49-(0)30-1208820-90
Email info@picoquant.com
Web www.picoquant.com