




C-Slide

*The New Name
for Countess Slides*




For use with Automated Cell Counters



C-Slide

Cell Counting Chamber Slides

- No need for coverslips
- Disposable, eliminates the need for washing
- Reduce exposure to infectious samples and risk of contamination
- Guaranteed reproducibility and reliability
- Made of sturdy quartz grade optical plastic



NanoEnTek

The developer & manufacturer
of Countess and Countess slide.

Specifications

Sample volume : 10 μ L

Staining method : Trypan blue stain

Dimension : 75 (L) x 25 (W) x 1.8 (H) mm



For use with Automated Cell Counters



Countess



Countess II

Ordering information

Catalog Number	Product Name
CS-050	C-Slide, Cell Counting Chamber Slides 50 slides (100 tests)

NanoEnTek

e-mail sales@nanoentek.com
website www.nanoentek.com

NanoEnTek, Inc.

851-14, Seohae-ro, Paltan-myeon, Hwaseong-si,
Gyeonggi-do, 18531, Korea
Tel : +82-2-6220-7940 / Fax : +82-2-6220-7999

NanoEnTek America, Inc.

240 Bear Hill Road, Suite 101, Waltham,
MA 02451, USA
Tel : +1-781-472-2558 / Fax : +1-781-790-5649

EC REP

MT Promedt Consulting GmbH
Altenhofstrasse 80,
66386 St. Ingbert, Germany